













SARSTEDT S-Monovette® System
TUBE GUIDE FOR GENERAL PRACTITIONERS



***TUBES MUST BE USED IN THE FOLLOWING ORDER**

FAILURE TO COMPLY MAY COMPROMISE PATIENT RESULT & PATIENT SAFETY

Tube / Bottle	Tube Type / Order*	Tests / Comments
 Green, 3ml	Sodium Citrate Coagulation 1st	PT, APTT, Fibrinogen  BLOOD TUBE MUST BE FILLED TO THE GREEN FILL LINE...
 White, 4.9ml	Serum 2nd STAND UPRIGHT ONCE DRAWN	Viral / Bacterial Serology
 Brown, 4.9ml	Serum Gel 3rd STAND UPRIGHT ONCE DRAWN	All Immunology One separate full tube for Immunology LDH, Lithium Drug Levels
 Orange, 4.9ml	Lithium Heparin Gel 4th	Routine Biochemistry Profiles Routine Endocrine Tests One full tube sufficient for Biochemistry and Endocrine Tests
 Grey, 2.7ml	EDTA 5th	FBC, IM Screen Only tube accepted is SARSTEDT 05.1167.901 (CS02146 MMUH Euro Route Logistics)
 Pink, 2.7ml	EDTA 6th	HbA1c, PTH, Red Cell Folate, NT-Pro BNP All of the above require a separate tube Only tube accepted is SARSTEDT 04.1917.001 (CS02139 MMUH Euro Route Logistics)
 Yellow, 2.7ml	Fluoride EDTA 7th	Glucose, Glucose Tolerance Test (GTT) Only tube accepted is SARSTEDT 04.1918.001 (CS0143 MMUH Euro Route Logistics)
 Purple, 3.5ml	ESR 8th	ESR  BLOOD TUBE MUST BE FILLED TO THE FILL LINE AT TUBE NECK ANGLE...

Refer to Laboratory for more specialised tests.
Please display this beside your venepuncture equipment in your clinical area

